PPS-110 Series



Precision Positioner Stage Reference Manual

(Open and Closed Loop Versions)



PPS-110 Precision Positioner Stage Reference Manual

Rev 2.03

MICRONIX USA, LLC

Tel: 949-480-0538 Fax: 949-480-0538

Email: <u>info@micronixusa.com</u> <u>http://micronixusa.com</u>

Table of Contents

1.	Intro	oduction	3
	1.1		3
	1.2	Recommended Controllers	4
	1.3	Technical Data	4
	1.4	Load Characteristics	4
2.	Mod 2.1		5 5
3.	Prep 3.1	paring to Install the PPS-110 Installation Preparation	6 6
	3.2	Package Contents	6
4.	Insta 4.1	alling the PPS-110 PPS-110 Installation	7 7
	4.1.	1 General Mounting	7
	4.1.	2 X-Y Mounting	9
5.	Con 5.1	0 • • • •	10 10
	5.1.	1 Open Loop Installation & Wiring Diagram	10
	5.1.	2 Closed Loop/Encoder Installation & Wiring Diagram	11
	5.1	Vacuum Environments	12
	5.1.	1 Handling and Preparation	12
	5.1.	Open loop Installation & Wiring Diagram	12
	5.1.	3 Closed Loop/Encoder Installation & Wiring Diagram	13
6.	Dim 6.1		14 14
	6.2	PPS-110 External Stepper Motor Version	15
7.	Stac 7.1		16 16
8.	Sup 8.1	,	17 17
	8.2	Units and Conventions	17
	8.3	Accessories	17
Α.	Арр		18
	A.1	• • •	18
	A.1.	1 Standard Atmospheric DB-9 Male Motor Connector	18



PPS-110 Precision Positioner Stage

1 5 1	10 Treeision Fosilioner stage	— Reference Manua
A.1.2	Stepper Motor Specifications	18
A.2 Li	near Motor	19
A.2.1	Standard Atmospheric DB-9 Male Motor Connector	19
A.2.2	Linear Motor Specifications	19
A.3 P	iezo Motor	20
A.3.1	Standard Atmospheric DB-9 Male Motor Connector	20
A.3.2	2 Phase Piezo Motor Wiring Connections	20
A.4 Li	mit Switches	21
A.5 O	pen Loop Vacuum Wiring Diagram	22
A.5.1	Straight Through 9-Pin Feed Through	23
A.6 L	Jsing the Digital Encoder Module	24
A.6.1	36BEncoder Module Pin-out	24
A.6.2	Operating and Electrical Specifications	25
A.6.3	Output Signals & Signal Termination for A quad B, Index and limits	25
A.6.4	Resolution	25
A.6.5	MII 6000 Digital Wiring Diagram	26
A.6.6	Straight Through 25-nin Feed Through	27




1. Introduction

1.1 Product Description


The PPS-110 is a high-precision, long travel linear stage. It is available with piezo motor, linear motor, or stepper motor with ballscrew drives. Steel recirculating ball bearings allow for smooth motion, and high guiding accuracy for loads up to 25 kg. Closed loop encoder resolution of 2 nm is achievable with a high precision digital encoder for piezo motor version.

Features:


- Travel ranges of 1000 mm
- Load capacity up to 25kg
- Closed loop encoder resolution of 2 nm (Piezo motor)
- Closed loop encoder resolution of 50 nm (Stepper motor)
- Closed loop encoder resolution of 10 nm (Linear motor)



PPS-110 300mm with External Stepper Motor option



PPS-110 100mm with Piezo Motor option



PPS-110 300mm with Linear Motor option

1.2 Recommended Controllers


The following controllers are available from MICRONIX USA:

- MMC-110 Piezo
- MMC-200 Stepper

1.3 Technical Data

Motor option	PM-003		SM-008	/SM-009	LM-004
Speed, max. [mm/sec]	5 (MMC-110)		10 (Higher spe	ed on request)	250 (depending on load)
Encoder option	Open Loop	Digital (RS-422)	Open Loop	Digital (RS-422)	Digital (RS-422)
Resolution, typical [µm]	0.001	0.002	0.1	0.05	0.01
Repeatability, bi-directional [µm]	n/a	± 0.1	± 4	± 0.2	± 0.1
Repeatability, uni-directional [µm]	n/a	0.1	0.4	0.2	0.1


1.4 Load Characteristics



Load Characteristics	F _x [N]	F _y [N]	F _z [N]	M _x [Nm]	M _y [Nm]	M _z [Nm]	k _{ax} [µrad/Nm]	k _{ay} [µrad/Nm]
PM-003	5	150	250	60	30	35	50	50
SM-008/SM-009	30	150	250	60	30	35	50	50
LM-004	150 peak 50 continuous	150	250	60	30	35	50	50

2. Model configurations

2.1 PPS-110 Order Numbers



Contact MICRONIX USA for custom applications and stacking configurations.

^{*}Linear motor is not available with open loop and high vacuum environment.

[†] Only available with Piezo and Stepper motor option or externally mounted UHV stepper motor.

3. Preparing to Install the PPS-110

3.1 Installation Preparation

When mounting the stage, it is important to consider the flatness of the mounting surface, as the stage will conform to the shape of that surface. A surface that is not flat can adversely affect the performance and structural integrity of the stage.

The stage is calibrated and guaranteed to be within specification at 20°C ±5°C, unless otherwise specified. The operating conditions of the stage are as follows:

- Mount to a clean and flat surface which is free of debris, burrs, and dings
- An indoor atmosphere free of corrosive gases, excessive dust, and condensation
- Temperature range of 0-40°C
- Relative humidity between 20-80%
- Locate away from water, heat, and electrical noise

3.2 Package Contents

If product is damaged or there are missing components, contact MICRONIX USA immediately. Do not discard product packaging in case of return shipment.

Package Should Contain:

- PPS-110 Linear Stage
- Reference Manual
- Any other previously agreed upon components such as a controller




4. Installing the PPS-110

Mounting patterns require M4 socket head cap screws or pan head machine screws for mounting, as well as M2.5 x 6mm dowel pins for precision alignment. Additional brackets and screws may be required for custom applications.

4.1 PPS-110 Installation


4.1.1 General Mounting

For general mounting configurations, mount the base to the mounting surface using the thru holes. Move the carriage to access base mounting pattern. (Please note, it is possible to move the carriage of the linear and piezo motor configurations manually without damaging the stage, however, for *stepper versions the motor must be driven by a controller to reposition the carriage.)



Metric M4 Top Mount Mounting Pattern

*For internal motor configurations: Mounting screws must be flush with base surface.




Travel	Α	В
lidvei	[mm]	[mm]
100mm	195	195
200mm	195	195
300mm	325	195
500mm	585	455
1000mm	1105	975


4.1.1.2 Metric Mounting: M4 Bottom Mount

Reauires:

2 x M2.5 x 6mm Dowel Pins 4 x M5 Socket Head Cap Screws



Metric M5 Bottom Mount Mounting Pattern




Travel	Α	В
navei	[mm]	[mm]
100mm	195	195
200mm	195	195
300mm	325	195
500mm	585	455
1000mm	1105	975

4.1.2 X-Y Mounting

For additional mounting configurations see Section 7: Stacking Configurations. Please note that it is possible to move the carriage of the linear and piezo motor configurations manually without damaging the stage, however, for stepper versions the motor must be driven by a controller to reposition the carriage.

Requires:

2 x M2.5 x 6mm Dowel Pins 4 x M4 Socket Head Cap Screws (Low Profile recommended)




5. Connecting the PPS-110

5.1 Atmospheric Environments

For controller information refer to the appropriate MMC controller manual. Wiring diagrams are consistent for all piezo, stepper, and linear motor assemblies.

5.1.1 Open Loop Installation & Wiring Diagram


Connecting the PPS-110 in an open loop configuration only requires that the D-sub 9 Pin Motor Cable be connected to a compatible controller. No other cables or components are required. Please note, open loop configurations are only available for <u>stepper or piezo</u> motor versions (See Appendix A.1 or A.3 for pinout configurations).




5.1.2 Closed Loop/Encoder Installation & Wiring Diagram

Using the PPS-110 stage with an encoder requires a closed loop compatible controller that recognizes the proper type of encoder feedback. Connect the stage as shown below.

5.1.2.1 MII 6000 Digital Encoder Wiring Diagram






5.1 Vacuum Environments

5.1.1 Handling and Preparation

When preparing the stage for vacuum environments, take the necessary precautions (such as wearing gloves, clean room, clothing, etc.) when handling the stage as to avoid any contaminants. Maximum Bake-out temperature is 100°C. MICRONIX USA can supply the stage with vacuum compatible connectors: 9-pin female PEEK connector for open loop, 25-pin female PEEK connector for closed loop with MII 6000 digital encoder.

5.1.2 Open loop Installation & Wiring Diagram

Connecting an open loop PPS-110 in a vacuum chamber requires the use of a feed through connector at the vacuum chamber wall. The vacuum compatible PPS-110 will be supplied with wiring for a straight through feed through, not a cross over gender changer. MICRONIX USA supplies test connectors that simulate the vacuum feed through to allow for functionality testing prior to installation in a vacuum chamber. Note: Linear motor versions are not available for vacuum environments. For details regarding the pin-out and feed through specifications see the Appendix section A.5.




5.1.3 Closed Loop/Encoder Installation & Wiring Diagram

Closed loop installation of the PPS-110 stage in vacuum environments requires an intermediate feed through connector at the vacuum chamber wall that can accommodate both the motor cable, and the encoder cable.


The vacuum compatible PPS-110 will be supplied with wiring for a straight through feed through, not a cross over gender changer. MICRONIX USA supplies test connectors that simulate the vacuum feed through to allow for functionality testing prior to installation in a vacuum chamber. For details regarding the pinout and feed through specifications see the Appendix sections A.6.5, A.6.6.

5.1.3.1 MII 6000 Digital Encoder Wiring Diagram




6. Dimensions


6.1 PPS-110 Internal Motor Versions (Linear and Piezo)



Linear Motor Version




Travel [mm]	L [mm]	A [mm]	B [mm]
100	240	195	-
200	340	195	-
300	440	195	325
500	640	455	585
1000	1140	455	585



Piezo Motor Version

6.2 PPS-110 External Stepper Motor Version



Travel [mm]	L [mm]	A [mm]	B [mm]
100	240	195	-
200	340	195	-
300	440	195	325
500	640	455	585
1000	1140	455	585




7. Stacking Configurations

7.1 Configuration Examples


- > Additional configurations available upon request
- > Note: Stacking compatibility for all motor configurations.
- > Positioning according to:



No Adapters




X-Y 300x300mm SM




X-Y 200x200mm SM



X-Y 100x100mm SM



X-Y 300x300mm PM



X-Y 200x200mm PM



X-Y 100x100mm PM

8. Supplementary Information

8.1 Maintenance

- The PPS-110 series of linear stages utilizes a maintenance free design. Do not modify the stage or perform any maintenance unless specifically instructed to do so by MICRONIX USA personal. If the stage is not performing up to the original specifications, please contact MICRONIX USA.
- The PPS-110 series of linear stages are precision mechanical device and should be handled with care. Do not drop or mishandle the stage.
- Do not touch the bearings, as this will contaminate the lubrication and jeopardize the longevity of the stage.
- Follow the *Installation Preparation* requirements and use proper cable management to ensure a clean and safe operating environment.

8.2 Units and Conventions

All measurements in this document are in the metric system of units.

Metric Unit	English Unit
1 millimeter	0.0394 inches
1 micron	0.0000394 inches
1 Newton	0.2248 lbs.
1 Newton-meter	8.85 in-lbs.

8.3 Accessories

430828	430613
Z Bracket PPS-60	Z Bracket PPS-60
Used to adapt the PPS-110 series	Used to adapt the PPS-110 series
to a perpendicular PPS-60 for XZ	to a perpendicular PPS-60 for XZ
mounting configurations.	mounting configurations.

A. Appendix

A.1 Stepper Motor

A.1.1 Standard Atmospheric DB-9 Male Motor Connector

		Wire Color		
Pin	Function	External Stepper SM-003		
1	Motor A+	Red		
2	Motor A-	Blue		
3	Motor B+	Green		
4	Motor B-	Black		
5	Limit Switch GND	Brown		
6	Limit Switch+	White		
7	Limit Switch-	Violet		
8	N/C	N/C		
9	N/C	N/C		
* White moto	* White motor wire may be changed to green or blue in the future			

A.1.2 Stepper Motor Specifications

External Stepper Motor

Motor Type	2 Phase Bipolar
Phase Current	1.3 A Max
Step Angle	1.8°
Steps	200
Coil-Resistance	1.9 Ohms
Coil-Inductance	1.7 mH
Pitch	1 mm/rev
Resolution/Full step	5 µm

[•] For vacuum prepared wiring, blue may be substituted for green/white, and black for red/white.

A.2 Linear Motor

A.2.1 Standard Atmospheric DB-9 Male Motor Connector

	Linear Motor					
Pin	Function	Wire Color				
1	Α	Red				
2	В	Black				
3	С	White*				
4	N/C	N/C				
5	Limit Switch GND	Brown				
6	Limit Switch+	White				
7	Limit Switch-	Violet				
8	N/C	N/C				
9	N/C	N/C				

A.2.2 Linear Motor Specifications

Linear Motor


Continuous Force	2.7 N
Continuous Current	0.8 Arms
Acceleration Force	10.8 N
Acceleration Current	3.4 Arms
Force Constant (K _f)	3.2 N/Arms
Back EMF (Ke)	1.1 V/m/s
Resistance 25°C	6.8 Ohms
Inductance	1.0 mH
Electric Time Constant	0.147ms
Rated Voltage (AC)	240 V
Fundamental Motor Constant (Km)	1.23 N√ <i>W</i>
Magnetic Pitch (North-North)	30mm

A.3 Piezo Motor

A.3.1 Standard Atmospheric DB-9 Male Motor Connector

Piezo Motor					
Pin	Function Wire Color				
1	Phase 1 + Motor 1	Re	ed		
1	Phase 1 + Motor 2	Wh	ite		
5	Phase 1 GND Motor 1	Blo	ick		
5	Phase 1 GND Motor 2	Grey			
2	Phase 2 + Motor 1	e 2 + Motor 1 White Gree			
2	Phase 2 + Motor 1	ase 2 + Motor 1 White Blue			
5	Phase 2 GND Motor 1 Green		een		
5	Phase 2 GND Motor 2	Blue			
6	Limit +	White	Violet		
7	Limit - Violet		olet		
8	Limit GND	Brown			

A.3.2 2 Phase Piezo Motor Wiring Connections



A.4 Limit Switches


The limit switches are normally closed (when not activated) and should be connected to a compatible controller that recognizes these settings. A hard stop is designed into the PPS-110 body which will prevent the carriage from running away should the limit switches fail.

The mechanical limit switches are factory calibrated to ensure advertised travel length and cannot be adjusted by the customer. Optical limit switches can be factory adjusted to custom, specific travel lengths. For custom travel lengths, please contact MICRONIX USA.

A.4.1.1 Mechanical Limit Switches

Contact Rating	100 mA @ 30 V		
Contact Type	Normally Closed		
Operating Temperature	-25 to +70 °C		

A.4.1.2 Limit Switch Schematic






A.5 Open Loop Vacuum Wiring Diagram

Standard Cable Descriptions:

- A. PPS-110 Motor Cable Vacuum Side (Female Dsub 9 Pin Peek Connector)
- B. Atmospheric Motor Connector (Female Dsub 9 Pin to Male Dsub 9 Pin)

Wiring Diagram:





Piezo Motor Connector Pinout

1 1020	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,					
Description	Color	L1	L2	L3	L4	L5
Phase 1	Red	5	5	1	1	1
Phase 2	Yellow	4	4	2	2	2
Limit Switch -	Violet	8	8	7	7	7
Motor Ground	Green&Blk	1	1	5	5	5
Limit Ground	Brown	7	7	8	8	8
Limit Switch+	White	9	9	6	6	6

Internal Stepper Motor Connector Pinout


Description	Color	L1	L2	L3	L4	L5
Motor Phase A+	Green	5	5	1	1 (Green)	1
Motor Phase A-	Green & White	4	4	2	2 (White - Green TP)	2
Motor Phase B+	Red & White	3	3	3	3 (Black)	3
Motor Phase B-	Red	2	2	4	4 (Red)	4
Limit Ground	Brown	1	1	5	5 (Brown)	5
Limit Switch +	Violet	9	9	6	6 (Violet)	6
Limit Switch-	White	8	8	7	7 (White - Violet TP)	7
Shield	-	6	6	9	9 (Shield)	Casing






Female Dsub9 Connector - Rear View

Male Dsub9 Connector - Rear View



A.5.1 Straight Through 9-Pin Feed Through



Over Index

Mark

Over Left Limit

Mark

A.6 Using the Digital Encoder Module

When using the diaital external encoder configuration, the Encoder Module should display two green LED's indicating a power source and proper encoder alignment. A Red or Yellow Signal Level indicates LFD misalignment of the Encoder Head, if this occurs contact MICRONIX USA. Do not adjust the Encoder Head or scale. For more information refer to MicroE Systems Mercury Encoders.

Optimal Signal

Strength -

Marginal Signal

Strength-

A.6.1 Encoder Module Pin-out


Pin	Description	Pin	Description
1	*Right Limit+	9	Ground
2	Ground	10	*Left Limit+
3	*Right Limit-	11	*Left Limit-
4	Index-	12	Index+
5	B-	13	B+
6	A-	14	A+
7	+5V	15	(not used)
8	+5V		

^{*-}Limits must be specified at the time of order, and calibrated at the factory.

Note: Tri-state alarm: A and B are tri-stated if the encoder signal becomes too low for reliable operation.

Normal Operation






lercury 11 6000



A.6.2 Operating and Electrical Specifications

Power Supply	5VDC ±5% @ 140mA (No outputs terminated) @ 180mA (A, B, I, and both limits terminated); 50mA at the sensor
Operating Temperature	0 to 70°C
Humidity	10 - 90% RH non-condensing

A.6.3 Output Signals & Signal Termination for A quad B, Index and limits



^{*}Output signals are differential. Inverse signals are not shown for clarity.

A.6.4 Resolution

All closed loop stages are supplied with 20µm scales. The digital encoder module interpolates to a higher resolution as specified in the order. With a digital encoder an MMC controller has an achievable resolution of 2nm.



^{**}Note: At some interpolations values the index pulse may be aligned with other states of A or B than the ones shown.


^{***}Above are with reference to the sensor's optical centerline

A.6.5 MII 6000 Digital Wiring Diagram

Standard Cable Descriptions:

- A. PPS-110 Motor Cable Vacuum Side (Female Dsub 25 Pin Peek Connector)
- B. PPS-110 Encoder Cable Vacuum Side
- C. Atmospheric Motor Cable (Female Dsub 25 Pin to Male Dsub 9 Pin)
- D. Atmospheric Encoder Module Cable (Female Dsub 25 Pin to MII 6000 Interpolator Module)
- E. Encoder Module Adapter Cable (Female Dsub 15 to Female Dsub 9 Pin)

Wiring Diagram:




7 (Shield)

◆ Note: For the pinout of cable E, refer to the appropriate MMC manual.

Shield





◆ Note: For the pinout of cable E, refer to the appropriate MMC manual.

A.6.6 Straight Through 25-pin Feed Through

